



Case Study - Frome Community Hospital

The New Frome hospital opened on 2 June 2008 was built to replace the old Frome Victoria hospital. Patients are admitted from district general hospitals and from home and receive their medical care from their own GP practice.

The Challenge

The new Frome hospital required full installation of Patient TV equipment to the newly refurbished wards.

From the outset clinical staff wanted to minimise wall mounted items and therefore identified ceiling mounts as the solution.



The Solution

The T2 Elite series of ceiling mounted hospital LCD TV/monitor mounting solutions is designed for ease of use by both patients and clinicians alike. Each T2 Elite ceiling mount provides 360° of screen swivel and 90° downward tilt. This enables patients to view TV whilst sitting or lying in bed and also whilst seated in an adjacent chair. As a fully cable managed monitor mounting solution, the arm is ideal in hospital and dental environments.

Designed to meet modern healthcare infection control standards and to be cleaned with hospital grade /chlorine based cleaning solutions, making them ideal for use in clinical areas including operating theatres, recovery, dental and wards.

ErgoMounts worked closely with the architect and main contractor Laing O'Rourke and carried out a full site survey and assessment of the new wards and it was agreed that the T2 Elite series would both benefit and exceed all the criteria.

ErgoMounts carried out full installation of the 28 solutions and an ongoing maintenance contract has been agreed. All was achieved on time and on budget.

The newly rebuilt and refurbished hospital was officially opened in June this year by the Countess of Wessex.

ErgoMounts Limited

Unit 10 Pegasus Court North Lane
Aldershot Hampshire GU26 6QP

T +44 (0)1252 333326 **F** +44 (0)1252 330006

www.ergomounts.co.uk